

TransIT-HeLaMONSTER® Transfection Kit

Quick Reference Protocol

Instructions for MIR 2900, 2904, 2905, 2906, 2910S

Full protocol, SDS and Certificate of Analysis available at mirusbio.com/2900



SPECIFICATIONS

Storage	Store both TransIT-HeLa Reagent and MONSTER Reagent tightly capped at -20°C. Before each use , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

► PLASMID DNA TRANSFECTION PROTOCOL



Full protocol and additional documentation available at mirusbio.com/2900

Fill in volumes below based on culture vessel used for transfection (Table 1).

A. Plate cells

1. Plate cells in ___ml complete growth medium (per well).
2. Culture overnight. Cells should be ≥80% confluent at the time of transfection.

B. Prepare TransIT-HeLa:MONSTER:DNA complexes

1. Warm TransIT-HeLa and MONSTER Reagents to room temperature and vortex gently.
2. Place ___μl of OptiMEM® I Reduced-Serum Medium in a sterile tube.
3. Add ___μl plasmid DNA. Mix gently by pipetting.
4. Add ___μl of TransIT-HeLa Reagent. Mix gently by pipetting.
5. Add ___μl of MONSTER Reagent. Mix gently by pipetting.
6. Incubate at room temperature for 15-30 minutes.

C. Distribute complexes to cells

1. Add TransIT-HeLa Reagent:MONSTER:DNA complex mixture drop-wise to different areas of the well.
2. Gently rock plate for even distribution of complexes.
3. Incubate 24-72 hours.
4. Harvest cells and assay as required.

Table 1. Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm ²	3.8 cm ²	9.6 cm ²
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 μl	100 μl	250 μl
DNA (1 μg/μl stock)	0.5 μl	1 μl	2.5 μl
TransIT-HeLa Reagent	1.5 μl	3 μl	7.5 μl
MONSTER Reagent	1 μl	2 μl	5 μl

► Transfection Optimization

Determine the best TransIT-HeLa:DNA and MONSTER:DNA ratio for each cell type. Start with 3 μl of TransIT-HeLa Reagent per 1 μg of DNA. Vary the amount of TransIT-HeLa Reagent from 2-4 μl per 1 μg DNA to find the optimal ratio. Vary the amount of MONSTER Reagent from 0-5 μl per 1 μg of DNA.

For additional optimization tips, see [full protocol](#).



Reagent Agent®

Reagent Agent® is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: [mirusbio.com/ra](https://www.mirusbio.com/ra)

©1996-2021 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.

For terms and conditions, visit www.mirusbio.com

Rev0 100621