


pLIVE™ Vector*
Part of Product # MIR 5420, 5520, 5620, and 5320



FEATURE	VECTOR LOCATION (bp)
Mouse Alpha Fetoprotein (AFP) Enhancer II	8 - 226
Mouse Minimal Albumin Promoter	269 - 586
Intron 1	602 - 900
Multiple Cloning Site (MCS) (<i>Nhe I/Asc I/Sal I/Sma I/BamH I/Sac I/Sac II/Xho I</i>)	908 - 961
Intron 2	988 - 1329
Polyadenylation Region	1472 - 1805
Kanamycin Resistance Gene	1888 - 2683
pUC origin (ori)	2746 - 3334

* Patent pending.