


pLIVE™-SEAP Vector*
Part of Product # MIR 5620 and 5320



FEATURE	VECTOR LOCATION (bp)
Mouse Alpha Fetoprotein (AFP) Enhancer II	8 - 226
Mouse Minimal Albumin Promoter	269 - 586
Intron 1	602 - 900
Human Placental Secreted Alkaline Phosphatase (SEAP) Open Reading Frame (ORF)	933 - 2489
Intron 2	2557 - 2898
Polyadenylation Region	3041 - 3374
Kanamycin Resistance Gene	3457 - 4251
pUC origin (ori)	4315 - 4903

* Patent pending.